Data Sheet



- COL G -Recombinant Collagenase class I



| Item No. | Item Description |
|----------|------------------|
| 001-001 | COL G, 75 U |
| 001-002 | COL G, 300 U |
| 001-003 | COL G, 750 U |

- COL H -

Recombinant Collagenase class II



| Item No. | Item Description |
|----------|------------------|
| 002-001 | COL H, 750 U |
| 002-002 | COL H, 3000 U |
| 002-003 | COL H, 7500 U |

1. DESCRIPTION

COL G and **COL H** are recombinant collagenases (metalloproteinases) class I and class II respectively [1]. COL G and COL H are synthesized separately from *C. Histolyticum* genes by DNA recombination in *E. Coli BL21 AI* strain, bearing a Maltose Binding Protein (MBP) tag at the N-terminal end [2].

COL G and COL H are affinity chromatography purified proteins, highly pure, highly stable, lot-tolot consistent, endotoxin-free (\leq 10 EU/mg, LAL assay) and animal-free.

| CAS: | 9001-12-1 |
|-------------|---------------------------------|
| EC: | 3.4.24.3 |
| Grade: | Research Premium Grade |
| Form: | Lyophilized white powder |
| Quality: | Amylose Affinity Chromatography |
| Inhibitors: | EDTA, EGTA, Cys, Hys, DTT, 2- |
| | mercaptoethanol |
| Activators: | Ca ²⁺ |

Their molecular weights are ~135 kDa (COL G) and ~158.5 kDa (COL H). COL G and COL H are soluble in water or aqueous buffers and express their maximum activity at **pH 8**.

2. SUBSTRATES

COL G and COL H play different synergic roles in collagen digestion. Indeed, COL G expresses a higher activity against **native collagen**, specifically hydrolyzing **3D-helix regions**, while COL H expresses a lower activity against the 3D helix and a higher activity against **linear collagen regions** at the motif Pro-Y-Gly-Pro [3,4]. The mix of COL G and COL H expresses a **synergic activity** that results in efficient collagen digestion [5].

For tissue dissociation, protease addition is needed to hydrolyze non-collagenous proteins and other macromolecules present in the extracellular matrix [6].

3. ENZYMATIC ACTIVITY

COL G \geq 3.0 Units/mg^{*} **COL H** \geq 30.0 Units/mg^{*}

*according to Grassmann, one Unit liberates 1 μmol of Gly-Pro-Ala from Carbobenzoxy-Gly-Pro-Gly-Gly-Pro-Ala-OH (Fluka 27673) in 1 min at pH 7.4, 37 °C [7].

4. APPLICATIONS

For research use only.

Due to their high purity and specificity, COL G and COL H are especially indicated for the isolation of primary cells from liver, pancreas, heart, cartilage and stem cells from adipose tissue and others.

In these applications we recommend using a combination of COL G and COL H in a specific activity ratio, or according to the relevant isolation protocol in order to obtain an optimal collagen digestion in cell isolation. For other applications or suggestions, contact <u>info@abielbiotech.com</u> or visit <u>www.abielbiotech.com</u>.



5. PREPARATION METHOD

We recommend reconstituting the lyophilized COL G and COL H enzymes in the tissuedissociation buffer by injecting the **buffer directly into the vial**. Do not exceed an enzyme concentration of 30 U/ml (COL G) or 300 U/ml (COL H) to avoid precipitates.

Keep the vial on ice and periodically shake until the enzyme is completely dissolved. Filter with $0.22 \ \mu m$ mesh for sterility.

Prepare a mix of COL G and COL H solutions in a specific activity ratio and dilute according to your protocol working solution concentration.

Add protease to the mix at 4 °C according to the specific application. Thermolysin, pronase or neutral protease/dispase can be normally used. **Protease must be added immediately before use** to avoid catalytic processes in the enzymatic blend. The amount of protease will define the aggressiveness of your enzyme mixture. For suggestions about your specific protocol and application please contact <u>info@abielbiotech.com</u> or visit <u>www.abielbiotech.com</u>.

6. STORAGE AND STABILITY

Lyophilized COL G and COL H are stable at -80 °C up to two years. We recommend splitting in aliquots the reconstituted solutions at need and storing them at -20°C up to one month or -80°C up to 6 months.

To use aliquots later on, they can be diluted in reconstitutive buffer or can be directly added into the enzyme working solution.

▲Warning: We recommend avoiding multiple freezethaw cycles and exposure to frequent temperature changes.

REFERENCES

- [1] Matsushita, O. et al. (1999) J. Bacteriol. 181(3): 923–933.
- [2] Salamone, M. et al. (2012) Chem. Eng. Trans. 27: 259-264.
- Philominathan S.T. et al. (2009) J. Biol. Chem. 284(16): 10868-10876;
- [4] Matsushita, O. Et al. (1994) J. Bacteriol.176: 149-156
- [5] Breite, A.G. et al. (2011) *Transplant Proc.* 43(9) : 3171-3175
- [6] Salamone, M. et al. (2014) Chem. Eng. Trans. 38: 247-252.
- [7] W. Grassmann, et al, (1960) Z. Physiol.Chemie 322:267

For suggestions about your specific protocol or application of COL G and COL H, contact us:

Monica Salamone, PhD R&D Manager m.salamone@abielbiotech.com Alessandra Lo Cicero, PhD Product Manager a.locicero@abielbiotech.com Tel. +39 091 6615633 Fax +39 091 661 5628

Abiel s.r.l.- Viale delle Scienze, ed. 16 - 90128 Palermo - P. IVA 02398200812 Tel.: +39 091 661 5633 – Fax: +39 091 661 5628 - <u>www.abielbiotech.com</u> - <u>info@abielbiotech.com</u>